




Decision Tree

Which surface should you choose?

Get expert advice: 450 659-1991 – sales@leika.ca



Sizes 36", 39", 42", and 48" available upon request.


Our therapeutic surfaces

Descriptive guide: features and clinical indications

Mensa

Therapeutic surface
without pump

STAGE 1-4/HIGH RISK



Positioning

Prevention and full treatment,
without pump


Main features

- 70% reduction in pressure (silicone cells)
- Adapts to variable shapes and different body segments (polymer carved in cubes)
- Promotes blood circulation
- Excellent comfort and reduction of shear forces
- Firm surface: facilitates care and transfers

Airmersus

Therapeutic surface
without pump

STAGE 1-3



Positioning

Cloud effect: enveloping

Main features

- Self-adjusting therapeutic mattress system
- Optimal pressure relief, without pump
- Ideal for palliative care, chronic pain, bone metastases, and cachexia

K2M/K2MP

Dynamic surface

STAGE 1-UNDETERMINED/
MAXIMUM PERFORMANCE



Positioning

High performance/Complex cases

Main features

- Microprocessor technology adjusts pressure in real-time
- Multiple modes: alternating, static, maximum inflation and pulsation
- 24-hour inflation autonomy (remains inflated)
- Border option (K2MP)
- Ventilated cover
- Automatic pressure adjustment in seated position
- Individually deflating cells

Serene Air

Dynamic surface

STAGE 1-4/BARIATRIC PROFILE



Positioning

Comfort and simplicity

Main features

- Turbine motor, ventilated cells, low air loss
- Exceptional performance for bariatric patients
- Available in Edge option

Domus Auto

Dynamic surface

STAGE 1-4/MODERATE RISK/
EASY TO USE



Positioning

Safety and efficiency

Main features

- Auto-adjusting pump
- Easy to use
- Reliable for moderate-risk profiles

Discover our products



Get expert advice:
450 659-1991
sales@leika.ca



Sizes 36", 39", 42", and 48" available upon request.