



GELFLEX FOAM®  
available in several sizes

4"

MAXIMUM WEIGHT

**350**  
pounds

WARRANTY

**5**  
YEARS  
on mattress

**3**  
YEARS  
on cover

## Leika Adapt 2 with SRT - 36" ROHO Modular System

### TECHNICAL SPECIFICATIONS

The ROHO Adapt 2 modular system, thanks to its SRT (Shear Reduction Technology™), provides an innovative way to reduce the effects negative pressure, friction and shear inside the mattress. The firm side edges and 4-inch bolster fixed in GelFlex® foam at the foot of the mattress holds the foam sections hybrid and ROHO in place without snaps. Several configurations are possible by positioning the different sections of hybrid foam and ROHO depending on the patient's needs. It should be noted that SRT technology, usually located in the center of the mattress, has been moved to the head of the mattress for the Adapt 2.

#### ROHO Technology

Recognized for 50 years, ROHO technology still respects its 4 fundamental principles: low resistance for immersion, thrust vertical, moves in all directions, low friction and shear. It does not contain DEHP.

#### Microclimate cover and handles

The removable Recovery 5 cover with zipper from 360 degrees gives clinicians access to the interior of the mattress. Presence handles on the side of the mattress.

#### Shear Reduction Technology (SRT)

The layers of foam, not being glued over the entire length, slide over each other when the head of the bed is activated, thus absorbing the shear.

