



Discover the KOMP mattress, the therapeutic mattress designed to treat stage 1 to 4 pressure wounds for high-risk patients. The ability to adjust the alternating cells from low or high variations provides performance possibilities while maintaining optimal comfort.

Main Features

- Pressure adjustment: gradation to patient weight
- Available modes: Alternating, static, maximum inflation, seat inflation
- Cell design: Full
- Visual and audible alarms in case of low pressure or failure
- CPR function: Rapid deflation in emergency situations
- Remains inflated up to 12 to 24 hours without power
- Technical cover: Waterproof, breathable, antibacterial, fire-retardant
- Heel relief plugs
- Maximum capacity: 360 lbs (160 kg)
- M = Under-mattress
- P = Foam borders

Description

CONTROL UNIT (COMPRESSOR)

- Alternating mode every other cell: Optimal pressure redistribution
- Static mode: When mattress is used in immersion
- Raised backrest mode (Seat Inflate): Manual adjustment
- Maxfirm mode: Facilitates care and transfers
- Patient weight adjustment: Manual adjustment: gradation to patient weight

MATTRESS

- Size: 36" x 80" x 6"
- Inflation autonomy: 12 to 24 hours
- Cells: full

TOP COVER

- Quilted and ventilated cover, which allows more efficient maceration management
- No visible seams, latex-free

- Fire-retardant, water-repellent, antistatic, antibacterial/antimicrobial
- Protected zipper closure
- Cover options: Detachable and stretchable polyurethane cover, impermeable to liquids and permeable to water vapour

FIRM SIDE BORDERS

- Foam side borders

Other Information

ACCESSORIES

- Transport bag included
- Heel relief plugs included

