



Discover the Pro-Care Optima, the therapeutic mattress designed to treat stage 1 to 4 pressure wounds for high-risk patients. Intuitive automatic detection technology instantly adjusts pressure according to the patient to ensure comfort and safety at all times. The pump is excessively quiet and very simple to use.

Main Features

- Pressure adjustment: 3 comfort levels
- Available modes: alternating, static, maximum inflation, seat inflation
- Cell-in-cell design with heel relief for added safety
- Visual and audible alarms in case of low pressure or failure
- CPR function: rapid deflation in emergency situations
- Remains inflated up to 8 hours without power (transport mode)
- Technical cover: waterproof, breathable, antibacterial, fire-retardant
- Heel adjustment valve
- Maximum capacity: 36 x 80 x 8: 550 lbs (250 kg)

Description

CONTROL UNIT (COMPRESSOR)

- Alternating mode every other cell: Optimal pressure redistribution
- Static mode: When mattress is used in immersion
- Raised backrest mode (Seat Inflate): Adjustment is automatic
- Maxfirm mode: Facilitates care and transfers
- Comfort adjustment: Manual adjustment (soft, medium, firm)

MATTRESS

- Available sizes: 36" x 80" x 8"
- Weight capacity: 550 lbs (250 kg)
- Inflation autonomy: 8 hours in transport mode
- Cells: Cell-in-cell format

TOP COVER

- Detachable and stretchable polyurethane cover, impermeable to liquids and permeable to water vapour
- No visible seams, latex-free
- Fire-retardant, water-repellent, antistatic, antibacterial/antimicrobial
- Protected zipper closure

FIRM SIDE BORDERS

- Not available

Other Information

ACCESSORIES

- Transport bag included
- Option: Heel relief plugs
- Option: Integrated battery (12-hour operating autonomy)

