



Recyclable negative pressure wound therapy device.  
7-15-30 days

## Main Features

- Rechargeable negative pressure wound therapy: 7-15-30 day option
- Possibility to purchase kit, dressing, adapter, canister and pouch
- Compact and portable system, disposable at end of treatment
- Continuous treatment: 80 mmHg or 125 mmHg
- Intermittent with advanced variable pressure: 80/30 mmHg OR 125/30 mmHg
- 70 ml canister easy to change by user or caregiver
- Comes with carrying case, strap and bed support
- Dimensions: 111 X 76 X 57 mm
- Weight: 0.23 kg
- Complete treatment options:
  - Continuous negative pressure therapy
  - Lateral tension to incision line
  - Negative pressure environment facilitating additional circulation to incision site
  - Formation of cellular collagen, elastin and extracellular matrix
  - Aids wound healing
  - Suitable for dehiscent wounds, traumatic, partial thickness burns, flaps and grafts

## Description

### FOAM AND DRESSING PROPERTIES

- Foam dressing ready to apply with silicone surface, repositionable
- Absorption: foam dressing 6X its weight plus a 70 ml reservoir
- Dressing can be left in place up to 7 days
- Integrated multilayer dressing
- Polyurethane film
- Super absorbent layer
- Foam layer
- Silicone gel adhesive layer

### SAFE USE

- Protects incision from external contamination
- Holds incision wound lips together
- Helps reduce traction forces across incision
- Helps reduce edema
- Decreases fluid collection, thus reducing exudates and hematomas
- Helps reduce risk of wound dehiscence
- Prevents entry of external pathogens into wound site
- Reduces potential for surgical site healthcare-associated infections
- Non-traumatic removal (thanks to silicone layer)



## USER-FRIENDLY OPERATIONS FOR CLINICIANS

- Usage period up to 7, 15 and 30 days depending on prescribed treatment
- Facilitates additional circulation to incision site
- Dressing designed to be left in place for extended period
- Direct contact with sutures, staples, sterile adhesive strips or glue
- Audible and visual alarms and indicator for lifespan, blockages and leaks

## PUMP LIFESPAN

- Disposable (recyclable), choice of 7, 15 or 30 days

