

BEST-IN-CLASS DYNAMIC SURFACE FOR ADULTS & KIDS

K2M

Maximum protection.

- Active alternating pressure for high-acuity care.
- Pediatric PICU · Oncology · Palliative.



35-1000 lb

Weight capacity

1.7x

Up to 1.7× higher risk
in low-weight patients.

Jia YJ et al., J Tissue Viability, 2024

< 2 hr

Pressure injury risk
on inadequate surfaces.

CHU de Québec

Why low-weight patients are different



Thinner skin. Less subcutaneous tissue.
PI develops faster.



Adult surfaces are **oversized and miscalibrated**
for pediatric weight.



PICU and oncology stays create **sustained**
high-risk exposure.



Alternating pressure, low air-loss cover,
timed cycles offload continuously.



Calibrated for pediatric and palliative,
in addition to adult weight ranges.



Purchase or rental.

K2M for adults, pediatric & palliative patients



LEIKA MEDICAL EQUIPMENT

Toll free - 24/7
1 877 659-1991

Direct
450 659-1991

Order / PO
sales@leika.ca



Scan here to rent or buy